

## Air Transport Certification Ergo 07 L

ErgoConcept certifies that the Ergo 07 L is suitable for air transport.

The company ErgoConcept has based on the international regulation of the transport of lithium batteries that prohibits the transport of lithium batteries whose energy exceeds 300 Wh or the Ergo 07L wheelchair's battery has an energy of 192Wh.

The watt hour energy of the wheelchair battery is 192 Wh calculated in the following manner from voltage (V) and amperage (Ah):

Energy (Wh) =  $2 V \times 0.5 Ah = 1 Wh$ 

The Ergo 07 L battery is easy to dismantle, its energy is 24V and 8Ah.

Dimensions in centimeters: 52x12x3

Weight: 1.4 Kg

Plastic wrap

## Characteristics of the Wheelchair:

Brand: ErgoConceptModel: Ergo 07L

Dimensions folded in centimeters: 70x56x29

- Weight: 22kg

- Material Chassis: Aluminium